

6. Accessories

Part Number	Description
904-0013	Glass Dish

7. Technical Specifications

This equipment is intended for indoor use and will meet its performance figures within the ambient temperature range of 20°C to 30°C, with maximum relative humidity of 80% (non-condensing). Installation Category II (transient voltages). Pollution Degree 2 in accordance with IEC 664. Suitable for operation at altitudes of up to 6500 feet.

Specifications:

Temperature Set Range:	(Ambient +5°C) to 65°C
Operating Temperature Range:	20°C to 30°C
Temperature accuracy:	± 1.0°C at 55°C
Stability:	± 0.2°C
Supply Voltage Range:	
145702 and 145951	115V ± 10%, 60 Hz
145702-2 and 145951-2	230V ± 10%, 50/60Hz
Power Rating:	324W
Heating Rate:	Ambient to 55°C within 45 minutes

8. Maintenance and Service

All Boekel laboratory products are designed to comply with IEC1010-1. No routine maintenance is required. There are no user serviceable parts in this product.

8.1 Cleaning

Disengage power cord prior to cleaning. The glass dish has molded handles to allow removal for easy cleaning. Paraffin build-up should be avoided around the top area of the Tissue Bath. Paraffin and liquids should also be cleaned as quickly as possible from the black chamber area to prevent glass dish heating problems. The outer casing may be cleaned with water and a damp cloth. Do not submerge or immerse the Tissue Bath in water. Before using any cleaning or decontamination method except those recommended by the manufacturer, users should check that the proposed method would not damage the equipment.

8.2 Replacement of Fuses

There are two supply fuses located in the fuse drawer. To replace the fuses:

- Disconnect the unit from the power supply.
- Remove the plug from the socket in the back of the unit.
- Pull back on the fuse drawer.
- Pull out the fuse drawer.
- Check and replace with the correct fuses if necessary. The fuses must be ¼" x 1 ¼" fast acting, rated 4A and 250V.
- Push the fuse drawer back in. Reconnect unit to the power supply.