Boekel Scientific Plasma Thawer
Model 301000
Operation Manual

Model Number: 301000

TABLE OF CONTENTS

1. SAFETY	3
2. PRODUCT INFORMATION	4
2.1 INTRODUCTION	4
2.2 CONTRAINDICATIONS	4
3. ASSEMBLY	4
3.1 UNPACKING	4
3.2 INSTALLATION	4
4. OPERATION	5
4.1 CONTROLS AND DISPLAYS	5
4.2 MOTOR AGITATION SWITCHES	6
4.3 CIRCUIT BREAKER	6
4.4 TEMPERATURE CONTROL	6
4.5 OVER TEMPERATURE ALARM	6
4.6 THAWING CYCLE	6
4.6.1 Getting Started	6
4.6.2 Temperature Validation	7
4.6.3 Adjusting the Cycle Times	7
4.6.4 Loading the Boekel Scientific Plasma Thawer	7
4.6.5 Starting a Cycle.....	8
4.6.6 Pausing a Cycle	8
4.6.7 Cycle Completion	8
4.6.8 Lowering the Basket Assemblies for storage	9
5. SPECIFICATIONS & OPERATING CONDITIONS	9
6. RECOMMENDED MAINTENANCE	9
6.1 WEEKLY	9
6.1.1 Draining and cleaning the tank.	9
6.1.2 Cleaning the outer case	10
6.2 QUARTERLY	10
6.2.1 Temperature Calibration	10
6.2.2 Over Temperature High Alarm Test	11
8. WARRANTY AND SERVICE	14
8.1 WARRANTY	14
8.2 SERVICE	14

Indications for Use of the Boekel Scientific Plasma Thawer

The Boekel Scientific Plasma Thawer is intended for use in blood banks and laboratories where it is required to thaw up to four units of Fresh Frozen Plasma (FFP) or Cryoprecipitate AHF, Pooled Cryoprecipitate AHF, PF24RT24, or Plasma Cryoprecipitate Reduced Within 24 Hours After Phlebotomy (PF24). This device is not intended for the thawing of, but not limited to, Red Blood Cells, Hematopoietic Progenitor Cells, or Umbilical Cord Blood.

1. Safety

The following symbols marked on the equipment mean:

Caution: Please read and understand all necessary installation and operating instructions prior to operating the system.

Always observe the following safety precautions:

- Use only as specified by the operating instructions or intrinsic protection may be impaired. After transport or storage in humid conditions, dry out the unit before connecting it to the supply voltage. During drying out intrinsic protection may be impaired.

- **RISK OF TRAPPING:** Do not attempt to add or remove bags from the baskets until the agitation function has completely stopped.

- Connect only to a power supply with a voltage corresponding to that on the serial number label.

- Connect only to a power supply that provides a safety ground terminal.
- Disconnect power before servicing, moving, cleaning, filling, or draining.
- Do not connect to a power supply or switch on before filling the bath
- Do not block or restrict ventilation slots.
- It is the user's responsibility to carry out appropriate decontamination if hazardous material is spilled on or inside the equipment.

- Before using any cleaning or decontamination method, users should verify that the method of cleaning will not damage the unit.
- If the temperature displayed is more than $\pm 1^{\circ}\text{C}$ off from the calibrated independent thermometer, the unit requires service and should not be used.

- Protective earth conductor, removal of the conductor could impair the safety afforded to the product.

- Electromagnetic interference could affect the operation of the system if it is used in the vicinity of devices that have not been evaluated to the relevant EMC standard/s.

2. Product Information

2.1 Introduction

The Boekel Scientific Plasma Thawer offers both controlled temperature and agitation to thaw various volumes of Fresh Frozen Plasma (FFP) or Cryoprecipitate AHF, Pooled Cryoprecipitate AHF, PF24RT24, or Plasma Cryoprecipitate Reduced. Two independently controlled baskets can accommodate up to two plasma bags per basket. The digital microprocessor control offers rapid temperature recovery and optimal temperature stability to reduce the time required to thaw multiple bags. The easy to program user interface allows multiple time options for thawing various sized bags. Upon completion of a timed cycle, the baskets will automatically rise from the bath. Advanced safety features will notify the user of unsafe conditions and raise the baskets from the bath upon alarm activation.

2.2 Contraindications

The Boekel Scientific Plasma Thawer is not to be used for the thawing of, but not limited to, Red Blood Cells, Hematopoietic Progenitor Cells, or Umbilical Cord Blood.

3. Assembly

3.1 Unpacking

- Remove packing materials carefully and retain for future shipment or storage of the unit.
- Inspect for damage. Report all shipping damage to the carrier immediately. Shipping damage is covered by the carrier and repair/replacement for shipping damages must be coordinated through the carrier.
- Complete and return the Warranty Registration Card.
- Package should contain:
 - (1) Boekel Scientific Plasma Thawer
 - (1) Power line cord
 - (1) Operating instructions
 - (1) Easy-Read® 0-60°C Thermometer
 - (1) Bath Drain Tube
 - (1) Box Disposable Plasma Overwrap Bags (250)
 - (1) Sample Bottle Bath Stabilizer

3.2 Installation

1. Place the Boekel Scientific Plasma Thawer on a flat and stable surface, making certain the bath has at least 10" of clearance above the bath to allow sufficient clearance of the baskets during operation.

- There are two water fill lines located on the rear surface of the bath wall. Fill the bath with tap or distilled* water only, to the mid-line level between the two water fill lines. The following diagram represents the bath fill line orientation:

___ ___ ___ ___ Maximum Water Level Fill Line

___ ___ ___ ___ Minimum Water Level Fill Line**


*Deionized water should not be used due to the risk of corrosion. Use only tap or distilled water.

**It is recommended that the bath be filled to the Minimum Water Level Fill line if thawing large volume, 500 ml bags.

- Install the thermometer in the thermometer grommet located between the two baskets on the back wall of the bath.
- Fit the power line cord into the IEC power socket on the rear of the unit and plug the power cord into a properly grounded outlet.

4. Operation

4.1 Controls and Displays

- Power Key (⏻) – Controls Power to the unit
- Left Basket Cycle Time Up (▲) and Cycle Time Down (▼) – Controls the cycle time of the Left basket.
- Left Basket Cycle Time Display - Displays the time remaining during a timed cycle.
- Left Basket Start cycle (⌚) – Starts the thawing cycle for the left basket.
- Left Basket Raise/Lower (⬆️) – Raises or Lowers the left Basket.
- Temperature Display- Displays the current temperature of the bath.
- Alarm Silence (🔇) – Silences the Audible Alarm
- Right Basket Cycle Time Display - Displays the time remaining during a timed cycle.
- Right Basket Cycle Time Up (▲) and Cycle Time Down (▼) - Controls the cycle time of the Right Basket
- Right Basket Raise/Lower (⬆️) - Raises or Lowers the Right Basket.
- Right Basket Start cycle (⌚) - Starts the thawing cycle for the Right Basket.

4.2 Motor Agitation Switches

Motor agitation switches are in the rear of the unit and are shielded to protect against accidental activation/deactivation. The top switch controls the agitation function of the right basket, and the bottom switch controls the agitation function of the left basket. When the switches are in the OFF position the bath is fully functional except for the basket agitation function.

4.3 Circuit Breaker

The unit is equipped with two circuit breakers on the rear of the unit. A breaker can be reset by pressing the circuit breaker button.

4.4 Temperature Control

The temperature controller is factory programmed and preset at 36.0°C fixed temperature.

4.5 Over Temperature Alarm

The factory preset over temperature alarm value is fixed at 37.0°C and is not adjustable. The over temperature alarm will sound if the water temperature in the bath exceeds the bath temperature set point of 36.0°C by 1.0°C. If the bath temperature exceeds 37.0°C, an audible over temperature alarm will sound, the baskets will rise from the bath, and the temperature display will flash between “ALr” and the current bath temperature. The Cycle time displays will flash “Er” indicating that the cycle was stopped due to an Over Temperature condition. During the investigation of the Over Temperature Condition, the Alarm Silence (🔇) button can be pressed to silence the audible alarm. Upon resolution of the alarm condition, the Alarm Displays can be cleared by turning the unit off and then turning the unit back on.

4.6 Thawing Cycle

4.6.1 Getting Started

To power the unit ON, press and release the power key once.

Please note: Upon power start up a display test is run to demonstrate that all display components are functioning properly.

To power OFF the unit press and hold the power key until the unit turns off. In preparation for running a thawing cycle, confirm that the water level is filled per section 3.2, power the unit ON, and allow the water temperature to stabilize to the temperature set point of 36.0°C± 0.3° C.

NOTE: Activating the basket agitation cycles will allow the bath temperature to stabilize more efficiently, see section 4.6.5, Starting a Cycle.

4.6.2 Temperature Validation

After the bath temperature has stabilized for approximately 30 minutes, verify that the temperature on the temperature display reads $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$. In the unlikely event that the temperature display does not read $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$ after stabilization, please contact Boekel Scientific for service.

Insert an independent calibrated thermometer into the bath. Read and record the temperature value on the independent calibrated thermometer and verify that the temperature of the bath is $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$. If the temperature reading of the bath is not within the temperature set point tolerance of $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$, then please refer to section 6.2.1, Temperature Calibration.

4.6.3 Adjusting the Cycle Times

The cycle times for the left and right baskets operate independently and allow the user to set different time requirements for each basket. The times can be set in one minute increments from 0-60 minutes.

The Up and Down keys are programmed with a speed-up algorithm to reach the desired value quickly. Adjusting the cycle time beyond 60 minutes will put the cycle timer into HOLD mode and “Ho” shall be displayed. The HOLD mode will allow a user to run a thawing cycle for an unspecified period until the user chooses to end the cycle manually. The cycle time set point set by the user shall be stored in the controller once a cycle is started and can be adjusted as necessary by using the corresponding left or right Basket Cycle Time Up or Down arrow keys prior to starting a new timed cycle. If the unit is powered OFF, or if power is lost and recovered, the cycle time will default zero and the user will need to set the desired thaw time. See Table 4.6.3 for some general thaw time guidelines.

NOTE: Before starting a thawing cycle always verify the thaw times detailed on Table 4.6.3 for appropriate volume and product being thawed.

Table 4.6.3 Thaw Times Based on Plasma Volumes

Type	Approximate Volume	Average Thaw Times (Minutes)*			
		1-Bag	2-Bags	3-Bags	4-Bags
Cryoprecipitate AHF	17-22 ml	4.0	4.0	5.0	5.0
FFP	250 ml	11.0	11.0	12.0	12.0
FFP	300 ml	14.0	15.0	16.0	16.0
FFP	500ml	28.0	34.0	34.0	34.0
*Thaw times listed are based on a thawing temperature of 36.0°C . All values are approximate and may vary slightly based on actual bag volume and storage conditions.					

4.6.4 Loading the Boekel Scientific Plasma Thawer

Press the Raise/Lower (↕) key for the Left and/or Right side to raise the basket. Insert a frozen plasma bag into a plasma overwrap bag to prevent contamination and place it into the basket assembly. To secure the overwrap bag to the basket, slip the slot in the bag over the overwrap hanging tab located at the top of the basket assembly.

4.6.5 Starting a Cycle

Press the Left Start cycle (⏻) button to begin the timed thawing cycle for the left basket assembly or the Right Start cycle (⏻) button to begin the timed thawing cycle for the right basket assembly. The corresponding basket assembly will lower into the bath and will begin the agitation of the basket. When a cycle is active, the decimal point shall blink in the lower right corner of the cycle time display and the timer will begin counting down.

4.6.6 Pausing a Cycle

Pressing the Left or Right Raise/Lower key while a cycle is in process will pause the cycle timer, de-energize the agitator and raise the corresponding basket. The display will flash and indicate the remaining time of the paused cycle. During a paused cycle, the following actions can be taken:

1. To resume the cycle, press the Raise/Lower (↕) key once to lower the baskets and resume the cycle.
2. To increase or decrease the remaining cycle time, press the Cycle Time Up (▲) or Cycle Time Down (▼) key. Changes to the time during a paused function will not change the original set cycle time.
3. To reset the Cycle if the sample has thawed, press, and HOLD the Start cycle (⏻) button for 3 seconds until an audible beep is heard and the cycle is reset to the original cycle setting.

4.6.7 Cycle Completion

Upon the completion of a cycle, the cycle timer will count down to zero, agitation will stop, and the basket will rise from the bath. An audible tone will beep five times and the Cycle time will reset to the original user set Cycle time.

To remove the thawed Plasma bag from the basket assembly, gently pull the plasma overwrap bag from the basket assembly and remove the thawed plasma bag from the overwrap bag.

NOTE: In the unlikely event that the basket assembly does not raise from the bath at the end of a cycle, the user can remove the thawed plasma from the bath by gently grasping the plasma overwrap bag with one hand while gently opening the basket clamp down mechanism with their other hand. With the basket clamp down mechanism open, the plasma overwrap bag containing the thawed plasma can be lifted and removed from the bath. This step can be repeated until all bags have been removed from the bath. After all bags have been removed from the bath, contact

Boekel Scientific for service.

Users can now reload the system or allow the system to remain idle until there is a need for additional thawing.

NOTE: Always verify the appropriate cycle time settings for the volume of product to be thawed prior to starting a new Thaw cycle.

4.6.8 Lowering the Basket Assemblies for storage

1. Remove FFP or Cryoprecipitate AHF bags from the unit.
2. Press the Left Basket Raise/Lower (↕) key to lower the Left Basket assembly. Press the Right Basket Raise/Lower (↕) to lower the right basket assembly.
3. When both basket assemblies have been lowered into the bath, the optional accessory Chamber cover, Model D1905565 can be used to reduce water evaporation.

5. Specifications & Operating Conditions

Boekel Scientific Plasma Thawer Specifications and Operating Conditions	
Maximum bag capacity (up to 500 ml bags)	4
Ambient Temperature Operation Range	18.0-30.0°C
Temperature Set Point	36.0°C (Factory set and fixed and is not adjustable)
Temperature Accuracy	+/- 0.3°C
Chamber Volume	4.5 – 5.0 gallons
Chamber Material	Stainless Steel
Baskets	(2) Stainless Steel
Chamber Drain Time	Approx. 5 minutes
Exterior Dimensions (w x h x d)	18.5" x 16.25" x 14.75"
Overall Exterior Dimensions (w x h x d)	21.75" x 23" x 15.5"
Chamber Dimensions (w x h x d)	15.5" x 11" x 7.5"
Electrical	120V, 50/60Hz 7A, 900W
Net Weight	58 lb.
Storage Temperature Range	20°C - 60°C
Humidity	25 – 90%
Altitude	0 – 2000 m above sea level
Atmospheric Pressure	87-108 kPa

6. Recommended Maintenance

6.1 Weekly

- 6.1.1 Draining and cleaning the tank.

1. Power off the unit by pressing and holding the Power Key.
2. Before draining, disconnect the power cord from the power supply socket.
3. Place the open end of the drain hose into sink or drain.
4. Place the drain hose assembly into the drain port located on the left side of the unit. As soon as the drain hose assembly is snapped into the drain port the bath will begin to drain.
5. When the bath is drained, be sure to remove the drain hose assembly by pressing the button on drain port and remove hose from drain port.
6. Remove the basket assemblies from the bath.
7. Clean the inside of the bath and basket assemblies with a soft cloth, mild detergent, and water.
8. Install the basket assemblies.
9. Refill the bath chamber as described in section 3.2.

6.1.2 Cleaning the outer case

1. The outer casing can be cleaned with a mild detergent and cloth dampened with water.

6.2 Quarterly

6.2.1 Temperature Calibration

In the event that the temperature reading of the bath is not within the temperature set point range of $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$ per section 4.6.2 Temperature Validation, the following Calibration procedure can be utilized to calibrate the bath.

NOTE: In the event that the temperature of the independent calibrated thermometer exceeds $\pm 1.0^{\circ}\text{C}$ from the Temperature set point of 36.0°C , the unit will require service and should NOT be used.

1. The bath is brought up to the operating temperature of 36.0°C by the controller and allowed to stabilize with the agitators on, for approximately 10 minutes.
2. The difference in temperature between the temperature display and a calibrated independent thermometer placed in the center of the bath is noted.
3. Press and hold the Power key to power down the unit.
4. Press and hold both the Left Cycle Up and Left Cycle Down keys while pressing the Power key to enter the Calibration Mode.
5. When the Calibration Mode is active, the Temperature display will flash between "CAL" and the temperature of the bath.
6. The Right Basket Cycle Time Up and Right Basket Cycle Time Down keys are used to adjust the temperature to agree with the Temperature noted earlier in Step 2 from the independent calibrated thermometer that was previously placed in the bath.
7. The user presses the power key to exit the Calibration Mode and save the value. Users can revert to normal operation by pressing the power key quickly or can power OFF the unit by pressing and holding the power key.
8. Upon power ON, verify that the temperature display matches the

temperature of the calibrated independent thermometer. If a discrepancy exists between the two readings, the calibration procedure should be repeated.

6.2.2 Over Temperature High Alarm Test

1. Remove FFP or Cryoprecipitate AHF bags from the unit.
2. Verify that the bath temperature has stabilized and that the display reads $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$.
3. Ensure baskets are in the raised position before continuing.
4. Power OFF the unit by pressing and holding the power key until the unit powers off.
5. From a Power OFF state, press and hold the Right Cycle Timer Up, Right Cycle Timer Down, and Right Start Cycle keys and press the Power Key to activate the Over Temperature High Alarm Test mode.
6. Upon initiation of the Over Temperature High Alarm Test mode, the baskets will lower approximately 0.5 inches and the display shall flash between “tST” and the temperature value of the bath.
7. The unit heater will heat the bath water and when the bath temperature exceeds 37.0°C , the alarm shall sound, and the baskets shall return to a fully raised position.
8. “PAS” will be displayed if the unit passed the Over Temperature High Alarm Test.
9. “FAL” will be displayed if the unit failed the Over Temperature High Alarm Test. If the Unit has failed, contact Boekel Scientific for Service.
10. After performing the Over Temperature High Alarm Test, power OFF the unit by pressing and holding the power key and allow the bath temperature to stabilize to $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$ before use.

7. Troubleshooting

Symptom	Possible Cause	Action
Power		
Unit does not Power On	• Power cord not connected to Electrical receptacle on Unit, and/or Wall outlet.	• Fit the power line cord into the IEC power socket on the rear of the unit and verify if plugged into a properly grounded outlet.
	• Power supply failure	• Verify that power source is active and regulated.
	• Circuit Breaker has been activated.	• Reset circuit breaker.
	• Faulty fuse on main Circuit Board.	• Contact Boekel Scientific for Service
	• Malfunctioning Power Key	• Contact Boekel Scientific for Service
Temperature Control		
Temperature display does not stabilize to $36.0^{\circ}\text{C} \pm 0.3^{\circ}\text{C}$	• Faulty Temperature Sensor	• Contact Boekel Scientific for Service
Temperature Overshoot/Undershoot	• Faulty Temperature Sensor	• Contact Boekel Scientific for Service
Temperature display stabilizes but differs from	• Unit is in need of Calibration.	• Refer to section 6.2.1 Temperature Calibration

temperature noted on independent calibrated thermometer.		
Over Temperature Alarm		
Over temperature alarm sounds during start-up.	<ul style="list-style-type: none"> • Tank not properly filled with water. • Tank Filled with water. 	<ul style="list-style-type: none"> • Verify that the tank is filled with water per Section 3.2 • Verify that the water used to fill the tank was not warmer than 36.0°C.
Over temperature alarm sounds during use.	<ul style="list-style-type: none"> • System Malfunction 	<ul style="list-style-type: none"> • Contact Boekel Scientific for Service
Over temperature alarm does not sound during a high alarm test.	<ul style="list-style-type: none"> • System Malfunction 	<ul style="list-style-type: none"> • Contact Boekel Scientific for Service
Controls and Displays		
Display Characters are not illuminated.	<ul style="list-style-type: none"> • Display Malfunction 	<ul style="list-style-type: none"> • Utilize external clock or timer to monitor thaw time. • Contact Boekel Scientific for Service.
Control Keys not responding.	<ul style="list-style-type: none"> • Circuit Board Malfunction 	<ul style="list-style-type: none"> • Contact Boekel Scientific for Service
Display Message Received	<ul style="list-style-type: none"> • See Table 7.1 	<ul style="list-style-type: none"> • See Table 7.1
Basket Assembly and Agitation		
Basket Clamps will not open when basket is raised from the bath.	<ul style="list-style-type: none"> • Basket assembly is not installed correctly on agitator arm. 	<ul style="list-style-type: none"> • Re-install basket assembly on agitator arm.
The basket will not fully rise from the tank. Audible alarm still Functions.	<ul style="list-style-type: none"> • An overhead obstruction prohibits the raising of the basket from the tank. • Motor Failure 	<ul style="list-style-type: none"> • Remove overhead obstruction above tank. • Utilize external clock or timer to monitor cycle time. • Contact Boekel Scientific for Service
The basket will not fully lower into bath upon initiation of cycle. Audible alarm still functions.	<ul style="list-style-type: none"> • An obstruction prohibits the lowering of the basket into the tank. • Basket assembly is not installed correctly on agitator arm. • Motor Failure 	<ul style="list-style-type: none"> • Remove obstruction from tank. • Re-install basket assembly on agitator arm. • Contact Boekel Scientific for Service
Basket assembly will not agitate upon initiation of cycle. Audible alarm still functions.	<ul style="list-style-type: none"> • Motor Agitation switch turned off. • Motor Failure 	<ul style="list-style-type: none"> • Verify that the motor Agitation switch is on • Contact Boekel Scientific for Service.

Table 7.1 Control Panel Message Values

Left Cycle Display	Temperature Display	Right Cycle Display	Action
Er	ALR ←Flashing→ Temp C°	Er	Contact Boekel Scientific for Service
Er	ALR ←Flashing→ tHi	Er	Contact Boekel Scientific for Service

Er	ALR ←Flashing→ tLo	Er	Contact Boekel Scientific for Service
Er	cLc	Er	Contact Boekel Scientific for Service
Er	Temp C°	Cycle Time	Check for basket jams, clear jams, reset power, and attempt basket movement. If error message persists, contact Boekel Scientific for Service.
Cycle Time	Temp C°	Er	Check for basket jams, clear jams, reset power, and attempt basket movement. If error message persists, contact Boekel Scientific for Service.
	FAL		Contact Boekel Scientific for Service
88	888	88	Contact Boekel Scientific for Service
	PAS		Alarm Test Mode Passed- No Action Required
	tSt		Over Temperature High Alarm Test in Progress- See Section 6.2.2
	rEv	x.x	Display of firmware revision level- No Action Required
88	888	88	Display Test on Startup- No Action Required
01	CAL ←Flashing→ Temp C°		Calibration Invoked-See section 6.2.1

8. Warranty and Service

8.1 Warranty

When used in the appropriate laboratory conditions and in accordance with these operating instructions, Boekel Scientific warrants this product to be free of defective parts, material, and workmanship for a period of two years from the date of shipment. The liability of Boekel Scientific for any defective equipment during the warranty period shall be limited to the repair of defective equipment or replacement thereof without charge for parts or labor.

8.2 Service

A Boekel Scientific Returned Material Authorization (RMA) number provided by Boekel Scientific is required before any Boekel products are returned for any reason. Contact Boekel Customer Service at 1-800-336-6929. A Decontamination Certificate must be completed, signed by the user, and returned to Boekel Scientific prior to receiving the RMA number. Please be sure to mark the outside of the returned goods package with this RMA number to ensure prompt handling.

Boekel Scientific
855 Pennsylvania Blvd.
Feasterville, PA 19053 USA
Phone: (215) 396-8200 or (800) 336-6929
Fax: (215) 396-8264
Email: Boekel-info@boekelsci.com

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