

# iTrite Electronic Motorized Titrator

## iTrite

The iDrop Digital Burette delivers accurate, precise and convenient bottle-top titration, as well as optimum operator safety. It combines a magnetic stirrer with the titration process, making it easier for the operator to achieve satisfactory results in the field of analytical chemistry, food industry, water analysis, etc.



### Features

- High quality motor providing excellent precision and accuracy
- Electronic control decreases repetitive strain injury
- The control panel makes the operating process without the manual disturbance
- Built-in magnetic stirrer

Volume range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10µL
Volume accuracy	R=0.2% CV=0.07%
Velocity	16 stages
Operating temperature range	10°C~30°C
Quality control	DIN EN ISO 8655
Control type	External control, including the stirrer and the burette program