

# myBLOCK™ HL

With Heated Lid



- Prevents evaporation & condensation
- Compact, fits in the palm of your hand
- Exchangeable blocks, for tubes 0.2 or 0.5ml
- Independent temperature control of lid & block
- Room temperature to 37°C in 2 minutes

## Technical Data

Temp. Range:	Ambient +5 to 100°C
Temp. Accuracy:	+/- 0.3°C
Temp. Increments:	0.1°C
Temp. Uniformity:	+/- 0.2°C
Dimensions:	4 x 6 x 4.5 in./ 10.5 x 15 x 11.5 cm
Weight:	2 lbs / 0.9 kg
Electrical:	100 - 240V, 50-60Hz*

The new MyBlock™ HL is the first and only mini heat block featuring a heated lid for protecting samples against unwanted condensation and evaporation. Both the block and lid temperature can be selected independently in increments of 0.1°C, making this MyBlock the ideal choice for the most temperature sensitive applications, including immunoassays, LAMP assays, enzyme reactions & other hyper-sensitive temperature procedures.

Two blocks are available, the 0.2ml block (item BSH100-02) accommodates up to forty 0.2ml PCR tubes or up to 5 eight position PCR strips. The 0.5ml block (item BSH100-05) holds a maximum of twenty-four 0.5ml tubes.



## Ordering Information

**BSH200-HL\***

**MyBlock HL Mini Digital Dry Bath  
with heated lid (Blocks sold separately)**

### Available Blocks:

**BSH100-02**

**Block, 40 x 0.2ml or 5 x PCR strips**

**BSH100-05**

**Block, 24 x 0.5ml**

**BSH100-15HL**

**Block, 15 x 1.5 ml for BSH200-HL**

\* Includes US plug, for EU plug, please add (-E)

**Benchmark**   
Scientific

**PO Box 709  
Edison, NJ 08817  
Phone: 908-769-5555  
Fax: 908-222-1864**

**Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)**