

SERVICE: RETURN AUTHORIZATION IS REQUIRED

1. Should service be required, contact your salesperson or Boekel Scientific Customer Service at 1-800-339-6929 x5.

SPECIFICATIONS

Size (W x D x H):	21 ¾" x 21" x 24"
Weight:	60 lbs.
Chamber Size:	6.4 cubic feet
Shelf Size:	19 1/8" x 18 7/8"
Ground Lug:	¼" spade connector
Gas Inlet Port:	¼" O.D., 5/32" I.D.
Maximum Temperature	150°C

Boekel Scientific
855 Pennsylvania Blvd.
Feasterville, PA 19053
Boekel-Info@BoekelSci.com
1-800-336-6929 x5
Fax: 215-396-8264