

The right pipette for your applications

Transferpette® S, adjustable volume

Items supplied: Transferpette® S, DE-M marking, performance certificate, shelf / rack mount and silicone grease



	Capacity, µl (color-coded)		Description	A* ≤ ± %	µl	CV* ≤ %	µl	Subdivision µl	Cat. No.
●	0.1	– 1	D-1 ①	2	0.02	1.2	0.012	0.001	705868
●	0.1	– 2.5	D-2.5 ②	1.4	0.035	0.7	0.018	0.002	705869
●	0.5	– 10	D-10	1	0.1	0.5	0.05	0.01	705870
●	2	– 20	D-20	0.8	0.16	0.4	0.08	0.02	705872
●	5	– 50	D-50	0.8	0.4	0.3	0.15	0.05	705873
●	10	– 100	D-100	0.6	0.6	0.2	0.2	0.1	705874
●	20	– 200	D-200	0.6	1.2	0.2	0.4	0.2	705878
●	100	– 1000	D-1000	0.6	6	0.2	2	1	705880
●	500	– 5000	D-5000 ③	0.6	30	0.2	10	5	705882
●	1000	– 10000	D-10000 ④	0.6	60	0.2	20	10	705884

Transferpette® S, fixed volume

Items supplied: Transferpette® S, DE-M marking, performance certificate, shelf / rack mount and silicone grease



	Capacity, µl (color-coded)	Description	A* ≤ ± %	µl	CV* ≤ %	µl	Cat. No.
●	10	F-10	1	0.1	0.5	0.05	705808
●	20	F-20 ⑤	0.8	0.16	0.4	0.08	705816
●	25	F-25	0.8	0.2	0.4	0.1	705820
●	50	F-50	0.8	0.4	0.4	0.2	705828
●	100	F-100	0.6	0.6	0.2	0.4	705838
●	200	F-200 ⑤	0.6	1.2	0.2	0.2	705844
●	500	F-500	0.6	3	0.2	0.2	705854
●	1000	F-1000	0.6	6	0.2	0.2	705862
●	2000	F-2000	0.8	16	0.3	0.3	705864



All product information
can be found at shop.brand.de

Pipette-Sets, Transferpette® S, adjustable volume



Pipette Package 2

Type	Description	Cat. No.
MICRO	3 Pipettes (D-2.5, D-10, D-100), 3 x TipBox, 3 shelf/rack mounts	705890
MID	3 Pipettes (D-20, D-200, D-1000), 3 x TipBox, 3 shelf/rack mounts	705891
MACRO	3 Pipettes (D-1000, D-5000, D-10000), 3 shelf/rack mounts	705892
STANDARD	3 Pipettes (D-10, D-100, D-1000), 3 x TipBox, 3 shelf/rack mounts	705893
PIPETTE PACKAGE 1	5 Pipettes (D-2.5, D-10, D-20, D-200, D-1000), 3 x TipBox, bench top rack	705896
PIPETTE PACKAGE 2	5 Pipettes (D-10, D-100, D-200, D-1000, D-5000), 4 x TipBox, bench top rack	705897

Transferpette® S-8

Items supplied: Transferpette® S-8, DE-M marking, performance certificate, 2 x TipBox, filled with pipette tips from BRAND, 1 shelf / rack mount, 1 reagent reservoir, 1 mounting tool, silicone grease and 1 set of sealing rings made of silicone/FKM.



	Capacity, µl (color-coded)	Description	A* ≤ ± %	CV* ≤ %	Subdivision µl	Cat. No.
●	0.5 – 10	M8-10	1.6	1.0	0.01	705900
●	5 – 50	M8-50	0.8	0.4	0.05	705906
●	10 – 100	M8-100	0.8	0.3	0.1	705908
●	20 – 200	M8-200	0.8	0.3	0.2	705910
●	30 – 300	M8-300	0.6	0.3	0.5	705912

Transferpette® S-12

Items supplied: Transferpette® S-12, DE-M marking, performance certificate, 2 x TipBox, filled with pipette tips from BRAND, 1 shelf / rack mount, 1 reagent reservoir, 1 mounting tool, silicone grease and 1 set of sealing rings made of silicone/FKM.



	Capacity, µl (color-coded)	Description	A* ≤ ± %	CV* ≤ %	Subdivision µl	Cat. No.
●	0.5 – 10	M12-10	1.6	1.0	0.01	705920
●	5 – 50	M12-50	0.8	0.4	0.05	705926
●	10 – 100	M12-100	0.8	0.3	0.1	705928
●	20 – 200	M12-200	0.8	0.3	0.2	705930
●	30 – 300	M12-300	0.6	0.3	0.5	705932

Notice: When ordering instruments with DAkKS calibration certificates, the prefix 'DAkKS' must be added to the order number, e.g., DAkKS 705920.

* Calibrated to deliver (TD, Ex). Error limits according to the nominal capacity (= maximum volume) indicated on the instrument, obtained with instrument and distilled water at equilibrium with ambient temperature at 20 °C, and with smooth operation.

A = Accuracy, CV = Coefficient of variation

- ① D-1 for use only with nano-cap™ tip from BRAND.
- ② D-2.5: when using BRAND nano-cap™ tips or tips from other manufacturers, calibration and possibly adjustment is necessary.
- ③ D-5000 for use only with 5 ml tip from BRAND and 5 ml tip from Thermo Fisher Scientific Inc. (Finnpipette®).
- ④ D-10000 for use with 10 ml tip from BRAND and appropriate tips from Gilson® and Eppendorf®.
- ⑤ For use only with 2-200 µl pipette tips.