

Product Name:	PRF072WWW/0A
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator &amp; Auto Defrost Freezer Combo  5.2 cu. ft. (Refrigerator) / 1.3 cu. ft. capacity (Freezer)  Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible  Microprocessor temperature controller, refrigerator and freezer  Digital temperature display, refrigerator and freezer  Audible and visual high and low temperature alarms, refrigerator and freezer  2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer  Remote alarm contacts, refrigerator and freezer  Adjustable operating temperature range: Refrigerator: 1°C – 10°C;  Freezer:-15°C – -25°C  3 adjustable shelves (refrigerator)/2 adjustable shelves (freezer), shelves have back rails to prevent product from being pushed against back wall  Footed bottom shelf refrigerator, and freezer  Keyed door locks  Magnetic door gasket for positive seal  Integrated door handles  Leveling legs  Stacking kit  Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior  HFC-free refrigerant (R600a), foam insulation and packaging  Two year parts and labor warranty, plus an additional three year compressor parts warranty  Exterior dimensions: 23 3/4"W x 24"D x 53 1/8"H  Auto defrost freezer / cycle defrost refrigerator  115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer  UL/C-UL listed  Shipping weight: 215 lbs.  *Please note this unit requires 2 outlets.</p>