

Product Name:	PRF092WWW/0
Description:	<p>Corepoint Scientific™ White Diamond Series Refrigerator & Freezer Combination 9 Cu. Ft. capacity (4.6 Ref./4.2 Freezer) Two exterior solid swing doors right hinged* Microprocessor temperature controller, refrigerator and freezer Digital temperature display with sample simulated probe, refrigerator and freezer Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C Audible and visual high and low temperature alarms, refrigerator and freezer Remote alarm contacts, refrigerator and freezer 2 adjustable/1 fixed shelves refrigerator/2 fixed shelves freezer with guard rails on back in Refrigerator to prevent product from being pushed against back wall Back wall evaporator cover, Refrigerator Forced draft circulation Refrigerator Keyed door locks Magnetic door gasket for positive seal Recessed door handle, freezer; heavy duty door handle, refrigerator Leveling legs Stacking kit LED interior lights Two probe access ports (3/8"); one, refrigerator/one, freezer White powder coated exterior HFC-free refrigerant (R600a), foam insulation and packaging Three year parts and labor warranty, plus an additional two year compressor parts warranty Exterior dimensions: 23 ¾"W x 24 ½"D x 66 ¾"H Manual defrost freezer / cycle defrost refrigerator 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer UL/C-UL listed Shipping weight: 269 lbs. In stock, with next day shipping</p> <p>*Please note this unit requires 2 outlets. *This product is available left hinged, product number PRF091WWWLH/0</p>