

Rotor overview



Rotor for 48 x 1.5/2 mL microcentrifuge tubes.

Can be stacked using rotor spacer 022822101 to double the capacity to 96 x 1.5/2 mL tubes. (F-45-70-11; Cat. No. 022822047)



Rotor for 8 x 15 mL conical tubes.

(F-45-8-17; Cat. No. 022822225)



Rotor for 6 x 50 mL conical tubes.






Also accommodates 6 x 15 mL conical tubes. (F-35-6-30; Cat. No. 022822228)



Swing-bucket rotor for MTP, PCR and Deepwell Plates up to 27 mm high.

(A-2-VC; Cat. No. 022822241)

Ordering information

Tube/Plate type	Rotor	Max. tube dimensions in mm (Ø x L)	Max. tube capacity in mL	Catalog No.
	F-45-72-8 ²	8 x 31	72 x 0.5 mL	022822080 ²
	F-45-70-11	11 x 41 11 x 47	70 x 1.5 mL 70 x 2.0 mL	022822047
	F-45-48-11 ²	11 x 41 11 x 47	48 x 1.5 mL 48 x 2.0 mL	022822004 ²
	F-40-36-12	12 x 32	36 x 1.5 mL	022822209
	F-45-36-15	15 x 45-48 ¹	36 x 2-4 mL	022822128
	F-45-16-20	20 x 42-55 ¹	16 x 6.5/10 mL	022822136
	F-40-18-19	19 x 66	18 x 10 mL	022822161
	F-45-12-31	31 x 55	12 x 20 mL	022822217
	F-35-8-24	24 x 86-105 ¹	8 x 25 mL	022822187
	F-45-24-12	12 x 67-100 ^{1,3}	24 x 5-7 mL	022822144
	F-50-8-16	16 x 105-120 ¹	8 x 8-12 mL	022822233
	F-50-8-18	18 x 105-128 ¹	8 x 10-18 mL	022822179
	F-45-8-17	17 x 116-123 ¹	8 x 15 mL conical tubes	022822225
	F-36-6-30	17 x 116-123 ¹ or 30 x 116-123	6 x 15 mL conical tubes or 6 x 50 mL conical tubes	022822228
	A-2-VC	128 x 86 x 27	2 x MTP or skirted PCR plates 2 x Deepwell plates (DWP) (max. height = 27 mm)	022822241
		Work tray, for use with semi-/unskirted PCR plates, PCR strips and 0.2 mL PCR tubes, (set of 10) for use with A-2-VC rotor		951010031
		CombiSlide® adapter, for use of glass slide, (set of 2) for use with A-2-VC rotor		022638963

¹ Min./max. tube length.

² Rotor is stackable for double capacity using spacer 022822101.

³ Tube length = 75 mm: max capacity 24 tubes
Tube length > 75 mm: max capacity 12 tubes