

Product Name:	FF201WWW/0MHC
Description:	<p>Corepoint Scientific™ White Diamond Series Natural Refrigerant Flammable Storage Freezer  20 Cu. Ft. Capacity  One solid swing door, right hinged  Microprocessor temperature controller  Digital temperature display with sample simulated probe  Audible and visual high and low temperature alarms  Remote alarm contacts  Adjustable operating temperature range: -15°C – -25°C  Six fixed shelves  Keyed door lock  Magnetic door gasket for positive seal  Leveling legs  Probe access port (3/8")  Powder coat exterior  High density urethane foam cabinet and door insulation  HFC-free refrigerant (R290), foam insulation and packaging  Three year parts and labor warranty  Exterior dimensions: 30 ¾ "W x 28 ½" D x 74 5/8"H  <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i>  Manual defrost  115V, 60 Hz, 1.96 Amps, ½ HP  UL/C-UL listed  Unit meets NFPA guidelines 45 and 70  Shipping weight: 240 lbs.</p>