

---

## Specifications

Gel plate size (w × h)	SE600: 18 × 16 cm SE660: 18 × 24 cm
------------------------	--

---

Gel size SE600: 14 or 16 × 16 cm  
16 cm with 1 cm wide spacers  
14 cm with 2 cm wide spacers

---

Gel size SE660: 14 or 16 × 24 cm  
16 cm with 1 cm wide spacers  
14 cm with 2 cm wide spacers

---

Maximum watt 50 W

---

Maximum volt 1000 V

---

Maximum ampere 500 mA

---

Maximum temperature 45 °C

---

Environmental operating conditions Indoor use: 4–40 °C  
Humidity up to 80%  
Altitude up to 2000 m  
Installation category: II  
Pollution degree: 2

---

Dimensions SE600  
(w × h × d) 32 × 29 × 14 cm  
(12.5 × 11.5 × 5.5 in.)  
SE660  
32 × 37 × 14 cm  
(12.5 × 14.5 × 5.5 in.)

---

Product certifications EN61010-1,  
UL61010A-1,  
CSA C22.2 1010.1,  
CE Certified

---

**This declaration of conformity is valid only when the instrument is:**

- used in laboratory locations,
- used as delivered from Hoefer, Inc. except for alterations described in the user manual, and
- connected to other CE-labeled instruments or products recommended or approved by Hoefer, Inc.