

Specifications

Gel plate size (w × h)	SE600: 18 × 16 cm SE660: 18 × 24 cm
------------------------	--

Gel size SE600:	14 or 16 × 16 cm 16 cm with 1 cm wide spacers 14 cm with 2 cm wide spacers
-----------------	--

Gel size SE660:	14 or 16 × 24 cm 16 cm with 1 cm wide spacers 14 cm with 2 cm wide spacers
-----------------	--

Maximum watt	50 W
--------------	------

Maximum volt	1000 V
--------------	--------

Maximum ampere	500 mA
----------------	--------

Maximum temperature	45 °C
---------------------	-------

Environmental operating conditions	Indoor use: 4–40 °C Humidity up to 80% Altitude up to 2000 m Installation category: II Pollution degree: 2
------------------------------------	--

Dimensions (w × h × d)	SE600 32 × 29 × 14 cm (12.5 × 11.5 × 5.5 in.) SE660 32 × 37 × 14 cm (12.5 × 14.5 × 5.5 in.)
------------------------	--

Product certifications	EN61010–1, UL61010A-1, CSA C22.2 1010.1, CE Certified
------------------------	--

This declaration of conformity is valid only when the instrument is:

- used in laboratory locations,
- used as delivered from Hoefer, Inc. except for alterations described in the user manual, and
- connected to other CE-labeled instruments or products recommended or approved by Hoefer, Inc.