


## Operational



### SUB6 Mini Horizontal

The SUB6 Mini Horizontal gel electrophoresis unit is ideal for routine preparatory and analytical electrophoresis techniques.

#### Technical Specifications

Unit dimensions (W × L × H)	13 × 24 × 6.5 cm
Gel dimensions (W × L)	6 × 7.5 cm
Buffer volume	325 ml
Maximum sample capacity	32
Combs	2
Comb thickness	1, 1.5 or 2 mm
Comb throughput	8 to 16 samples
Comb slots	2
Migration distance between comb slots	3.5 cm
Recommended running voltage	70 to 90 V

### SUB6 Combs

Catalog No.	Thickness (mm)	Sample Throughput	Tooth Width (mm)	Max Spacing (mm)	Sample Volume in a 5-mm deep well (µl)
SUB6-C1-8	1	8	4.5	2	20
SUB6-C1-12MC	1	12*	2.5	2	11
SUB6-C1-16	1	16	2.2	1	10
SUB6-C1.5-8	1.5	8	4.5	2	30
SUB6-C1.5-12MC	1.5	12*	2.5	2	17
SUB6-C1.5-16	1.5	16	2.2	1	15
SUB6-C2-8	2	8	4.5	2	40
SUB6-C2-12MC	2	12*	2.5	2	22
SUB6-C2-16	2	16	2.2	1	20

\*Multi-channel pipette compatible