



designed for scientists



IKA PETTE MacroKit

/// Data Sheet

The IKA PETTE MacroKit pipette package covers a volume range of 100 - 10,000 µl, ideal for laboratories that are working in the milliliter volume range.

The package consists of three single channel pipettes with variable volumes in the macro-volume range, a complimentary stand that holds up to three pipettes and three pipette tip boxes. This makes it a convenient starter package for researchers.

FEATURES:

- Mechanical air-cushion pipettes manufactured in Germany



designed for scientists

- Fully automated calibration by robots
- Volume-locking feature to prevent any accidental volume changes
- Large and clearly visible display
- Ergonomic and lightweight pipette with low pipetting and tip ejection forces to reduce strain and fatigue
- Durable and robust construction using high quality UV resistant and chemical resistant materials
- Fully autoclavable to ensure sterility and compliance to stringent laboratory standards
- Maximum user convenience with tool-less disassembly and calibration adjustment for regular maintenance



designed for scientists

Technical Data

Design	Air displacement piston-operated pipette
Handling	mechanical
Volume	variable
Number of channels	1
Permissible ambient temperature [°C]	5 - 40