



Overview

Traps are connected between the evaporator and vacuum pump to provide additional protection to the pump.

Acid Trap Insert, contains acid neutralizing desiccant that changes from off-white to purple when exhausted. Chemical trap protects to the rotary vane pump from acid fumes.

Note that this item number includes the insert **ONLY**. Insert requires a Clear Polypropylene Canister with mounting bracket or Clear Polypropylene Canister with stand for proper use.

Catalog Number: 7814800


Specifications

- **Weight:** 4.0 lbs
- **Weight metric:** 1.8 kg
- **Required Accessories:** Chemical Trap Canister
- **Applicable Products:** CentriVap Concentrators, Precise Glove Boxes, RapidVap Vacuum Evaporation Systems

For Product Assistance



[KeiShaun Johnson](#)
Customer Coordinator



[Ryan Caton](#)
Sales Representative

