

The figure above: 1. Observation window; 2. Top cover; 3. Operation area; 4. Housing assembly; 5. Electrical control panel; 6. Rubber pad; 7. Hinge; 8. Rotor system; 9. Motor assembly; 10. Shock absorption system

1.1.3 Symbol and meaning

Symbol	Meaning
~	AC power
	Power on
O	Power off

1.2 Main parameter

This product belongs to class I medical device, and its basic performance parameters are as follows:

Main parameter	Specification
Power	AC 220V 50Hz
Power supply	150W
Max. speed	12000r/min
Speed control accuracy	±20rpm
Max. RCF	14170xg
Max. capacity	24 Capillary rotor
Timing range	1s-99min
Timing control accuracy	≤±1%
Noise	≤65 dB
Ambient temperature	5°C~40°C
Host machine size (L×W×H) mm	320*380*270 mm
Net weight	19kg
Gross weight	23kg

Table 1-1 Main parameter