











Rotors	C	SC	Cat. no.	Related centrifuge	Max speed (rpm)	Max RCF (xg)
				Thermo Scientific		
Benchtop rotors						
 Fiberlite F14-6x250 LE			75003662	Sorvall™ Legend™ XT, Multifuge™ X3, SL 40F series	10,000 11,000 ^[9]	15,317 18,533 ^[9]
 Fiberlite F15-6x100y			75003698	Sorvall Legend X1/XT, Multifuge X1/X3, SL 40F series	15,000	24,652
				Sorvall ST 16/40, Megafuge™ 16/40, SL 16/40 series	13,000	18,516
 Fiberlite F13-14x50cy	■		75003661	Sorvall Legend X1, Multifuge X1 series	8,500	12,359
			75006526	Sorvall Legend XT, Multifuge X3, SL 40F series	9,250 10,000 ^[9]	14,636 17,105 ^[9]
 Fiberlite F15-8x50cy	■			75003663	Sorvall Legend XT, Multifuge X3, SL 40F series	14,500
			Sorvall Legend X1, Multifuge X1		14,000 14,500 ^[9]	22,789 24,446 ^[9]
			75006516	Sorvall Legend T, Multifuge 3 series	12,000 14,500 ^[9]	16,741 24,446 ^[9]
 Fiberlite F21-48x2			75003664	Sorvall Legend X1/XT, Sorvall ST 16/40, Multifuge X1/X3, Megafuge 16/40, SL 16/40, SL 40F series	15,200	25,055
			75006527	Sorvall Legend T, Multifuge 3 series	15,000	24,400
 Fiberlite H3-LV			75003665	Sorvall Legend X1/XT, Sorvall ST 16/40, Multifuge X1/X3, Megafuge 16/40, SL 16/40, SL 40F series	3,600	2,738

9 With 230 V centrifuge. C = Conical tubes SC = Sample containment  Biocontainment certification by Public Health England, Porton Down, UK.

Thermo Scientific Fiberlite rotor adapters and accessories

Sorvall LYNX Superspeed rotors with Auto-Lock	Legacy Superspeed RC series rotors	Benchtop rotors	Ultraspeed rotors
Fiberlite rotor (Cat. No.)	Adapter description ⁽⁷⁾	No. of vessels per adapter	Cat. No.
39 mL Ultraspeed			
F50L-8x39 (096-087051)	13.5 mL Tube	1	010-1142
50 mL			
	30 mL Oak Ridge Tube	1	010-0167
	16 mL Oak Ridge Tube	1	010-0382
F20-12x50 LEX (096-124375)	15 mL Conical Tube	1	010-1123
F21-8x50y (096-084275)	10 mL Oak Ridge Tube	1	010-1306
F21-8x50y (46923)	10 mL BD Vacutainer™ Tube	1	010-1068
	3 mL BD Vacutainer Tube	1	010-1128
	1 mL BD Microtainer™ Tube	3	010-1127
50 mL Conical			
	50 mL Oak Ridge Tube	1	010-0377
F14-14x50cy (096-145075)	30 mL Oak Ridge Tube	1	010-1147
F13-14x50cy (46922)	16 mL Oak Ridge Tube	1	010-0376
F13-14x50cy (75003661)	15 mL Conical Tube	1	75100378
F13-14x50cy (75006526)	15 mL Millipore™ Filtration Device	1	010-1340
	10 mL Oak Ridge Tube	1	010-1311
	10 mL BD Vacutainer Tube	1	010-1124
100 mL			
	50 mL Oak Ridge Tube	1	75003102 ⁽¹⁰⁾
	50 mL Conical Tube	1	75003103 ⁽¹⁰⁾
	30 mL Oak Ridge Tube	1	75003094 ⁽¹⁰⁾
	16 mL Oak Ridge Tube	1	76002906 ⁽¹⁰⁾
F20-6x100 (096-064025)	15 mL Conical Tube	1	75003095 ⁽¹⁰⁾
F15-6x100y (75003698)	10 mL Oak Ridge Tube	1	75003093 ⁽¹⁰⁾
	10 mL BD Vacutainer Tube	1	010-1274
	6.5 mL Round Bottom Tube	2	750030921 ⁽¹⁰⁾
	3 mL BD Vacutainer Tube	3	010-1126
	1.5/2 mL Microtube	4	75003091 ⁽¹⁰⁾
	1 mL BD Microtainer Tube	6	010-1125
100 mL Ultraspeed			
F37L-8x100 (096-08056)	39 mL Tube	1	010-0189
	13.5 mL Tube	1	010-0191
250 mL			
	100 mL Oak Ridge Tube	1	010-1119
	50 mL Conical Tube	1	75100136
	50 mL Oak Ridge Tube	1	010-0138
F14-6x250y (096-062075)	30 mL Oak Ridge Tube	2	010-1072
F14-6x250y (78500)	16 mL Oak Ridge Tube	5	010-1074
F14-6x250 LE (75003662)	15 mL Corning™ Conical	5	75101073
	15 mL Conical Tube	5	010-1410
	10 mL Oak Ridge Tube	7	010-1309
	10 mL BD Vacutainer Tube	7	010-1117
	3 mL BD Vacutainer Tube	10	010-1138
500 mL			
	250 mL Conical Tube	1	010-1135
	250 mL Oak Ridge Tube	1	010-0151
	175 mL Nalgene Conical Bottle	1	010-0152
	100 mL Oak Ridge Tube	1	010-1114
	50 mL Conical Tube	1	010-1102
F12-6x500 LEX (096-062375)	50 mL Oak Ridge Tube	2	010-1112
F12-6x500 LEX (096-062185)	30 mL Oak Ridge Tube	3	010-1115
	16 mL Oak Ridge Tube	7	010-1105
	15 mL Conical Tube	6	010-1099
	10 mL Oak Ridge Tube	7	010-1308
	10 mL BD Vacutainer Tube	7	010-1103
	3 mL BD Vacutainer Tube	14	010-1137
1000 mL			
	500 mL Oak Ridge Tube	1	010-0145
	250 mL Conical Tube	1	010-1096
	250 mL Oak Ridge Tube	1	010-0150
	175 mL Nalgene Conical Bottle	1	010-1132
	100 mL Oak Ridge Tube	3	010-1093
	50 mL Conical Tube	5	010-0180
F9-6x1000 LEX (096-061075)	50 mL Oak Ridge Tube	7	010-0191
F10-4x1000 LEX (096-041075)	30 mL Oak Ridge Tube	7	010-1095
F10-4x1000 LEX (096-041053)	16 mL Oak Ridge Tube	15	010-1087
	15 mL Conical Tube	12	010-1079
	10 mL Oak Ridge Tube	18	010-1307
	10 mL BD Vacutainer Tube	18	010-1415
	6 mL BD Vacutainer Tube	22	010-1416
	4 mL BD Vacutainer Tube	19	010-1418
	2 mL Filtration Tube and 1.5 mL Conical Tube	12	010-1417
	1.8-2.7 mL BD Vacutainer Tube	30	010-1419
H3-LV Rotor			
	Promega™ Slicprep™ 96 Device (4 per run)	2	018-029032
H3-LV (75003665)	Standard Microplates (28 per run)	14	018-029031
	2 mL Deep-well Microplates (8 per run)	4	018-029031



Adapters sold in sets of 2, unless otherwise indicated.
 7 Actual fill volumes may vary from nominal volume
 10 Sold in pack of one.